PHILIPS

HearLink

Teknik veriler

HearLink 1730 | 1530 miniBTE T

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Technical data

HearLink 1730 | 1530 miniBTE T

HearLink miniBTE T is a behind-the-ear hearing instrument of the Philips HearLink family, suitable for slight to severe hearing losses. Powered by AI sound technology, the HearLink miniBTE T includes our audiological features in SoundMap 2. Thanks to Bluetooth® Low Energy technology,

it directly connects to iPhone, iPad, Mac and Android™ devices. The miniBTE T comes with the miniFit thin tube system, which includes a wide variety of domes and custom moulds.

Hook

MNB T

Thin tube 1.3 mm

MNB T

Thin tube 0.9 mm

MNB T

Technical features

- Direct audio streaming¹
- Hands-free communication²
- Bluetooth Low Energy technology
- NFMI (Near-Field Magnetic Induction)
- Single push-button
- Telecoil
- miniFit thin tube
- Hydrophobic coating
- LED visual indicator

Accessories

- Philips HearLink 2 app
- Philips Remote Control
- Philips TV Adapter
- Philips AudioClip

General features

- Digital Programmable
- Automatic or Manual VolumeControl
- Maximum Output Control System
- MPO-Maximum Power Output
- GC-Gain Control
- AGC-Automatic Gain Control
- Noise Reduction
- Feedback Management
- Dual Microphone
- FM Compatible(with Telecoil)
- 4 Programs
- 48 Channels

For information on compatibility, please visit hearingsolutions.philips.com/compatibility

Operating and charging conditions

Temperature: +1°C to +40°C (34°F to 104°F)
Humidity: 5% to 93% relative humidity, non-condensing
Atmospheric pressure: 700 hPa to 1060 hPa

Transportation and storage conditions

Temperature and humidity shall not exceed the mentioned limits for extended periods during transportation and storage:

Transport

Temperature: -25°C to +60°C (-13°F to 140°F)
Humidity: 5% to 93% relative humidity, non-condensing
Atmospheric pressure: 700 hPa to 1060 hPa

Storage

Temperature: -25°C to +60°C (-13°F to 140°F)
Humidity: 5% to 93% relative humidity, non-condensing
Atmospheric pressure: 700 hPa to 1060 hPa

1) Hands-free communication is available on select devices

2) From iPhone, iPad, Mac and select Android devices

WARNING: No modification of this equipment is allowed.

Apple, the Apple logo, iPhone, iPad, Mac and the Mac logo are trademarks of Apple Inc. registered in the U.S. and other countries. Use of the Made for iPhone badge means that an accessory has been design to connect specifically to the Apple product(s) identified in the badge and has been certified by the developer to meet Apple performance standards. Apple is not responsible for operation of this device or its compliance with safety and regulatory standards. Android™ is a trademark of Google LLC.

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


Feature overview

	HearLink 1730	HearLink 1530
SoundMap 2		
Amplification		
Frequency Bandwidth	8 kHz	8 kHz
Low Frequency Enhancement	•	•
Frequency Lowering	•	•
Noise Control		
Transition	•	•
Directionality		
Pinna	•	•
Fixed Directional	•	•
Adaptive Directionality	•	•
Dynamic Directionality	•	—
AI Noise Reduction		
Noise Reduction Mode	3 options	2 options
Special noise management		
Soft Noise Management	•	•
Wind Noise Management	•	•
Transient Noise Reduction	2 options	—
Feedback canceller		
Strength control	•	•
SoundTie 2		
Direct streaming ¹	•	•
Hands-free communication ²	•	•
Binaural coordination		
NFMI	•	•
Binaural Volume and Program Change	•	•
Non-Telephone Ear Control	•	•
Programming options		
General	•	•
Fitting Bands	14	12
Environments	10	8
Manual Listening Programs	4	4
HiFi Music	•	—
Data Logging	•	•
Adaptation Manager	•	•
CROS compatibility	•	•

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1) From iPhone, iPad, Mac and select Android devices
 2) Hands-free communication is available on select devices

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

Technical information
Omnidirectional mode is used
unless otherwise stated.

Hook / Thin tube 1.3
— Acoustic input: 60 dB SPL
- - - Magnetic input: 31.6 mA/m

Thin tube 0.9
— Acoustic input: 60 dB SPL
- - - Magnetic input: 31.6 mA/m

	Hook	Thin tube 1.3 mm	Thin tube 0.9 mm
OSPL90, Peak (dB SPL)	132	128	123
OSPL90, 1600 Hz (dB SPL)	127	123	116
OSPL90, HFA (dB SPL)	126	122	118
Full-on gain, Peak (dB) ¹	63	59	59
Full-on gain, 1600 Hz (dB) ¹	55	56	52
Full-on gain, HFA (dB) ¹	55	55	52
Reference test gain (dB)	48	47	41
Frequency range (Hz)	100-7500	100-7500	100-7500
Telecoil output, 1 mA/m field (1600 Hz) (dB SPL)	86	88	87
Telecoil output, 10 mA/m field (1600 Hz) (dB SPL)	106	*	*
Total harmonic distortion (Input 70 dB SPL, 500 Hz) (%)	<4	<5	<3
Total harmonic distortion (Input 70 dB SPL, 800 Hz) (%)	<4	<2	<2
Total harmonic distortion (Input 70 dB SPL, 1600 Hz) (%)	<2	<2	<3
Equivalent input noise level, Omni (dB SPL)	18	15	19
Equivalent input noise level, Dir (dB SPL)	28	*	*
Battery consumption, Typical (mA) ²	1.9	2.0	2.0
Battery consumption, Quiescent (mA) ²	1.9	1.9	1.9
Battery life, artificial measurement, hours ³	95	90	90
Expected battery life, hours (battery size 312 – IEC PR41) ⁴	50-55	*	*

* No measurement performed

- 1) Measured with the gain control of the hearing aids set to their full-on position minus 20 dB and with an input SPL of 70 dB. This is to obtain a gain response equal to the full-on gain response from e.g. IEC 60118-0:1983+A1:1994 but without influence of feedback.
- 2) Battery current is measured according to IEC 60118-0:1983/AMD1:1994 §7.11, IEC 60118-0:2015 §7.7 and ANSI S3.22:2014 §6.13 after a settling time of minimum 3 minutes.
- 3) Based on the standardised battery consumption measurement (IEC 60118-0:1983/AMD1:1994). The actual battery life depends on battery quality, use pattern, active feature set, hearing loss and sound environment.
- 4) Real usage battery life is shown as an estimated interval based on mixed use cases with variable amplification settings and variable input levels, incl. direct stereo streaming from a TV (25% of the time) and streaming from a mobile phone (6% of the time).

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HearLink 1730 | 1530 miniBTE T

2CC Coupler



Measured according to ANSI S3.22-2014,
IEC 60118-0:2015 and IEC 60318-5:2006

Technical information
Omnidirectional mode is used
unless otherwise stated.

Hook / Thin tube 1.3
— Acoustic input: 60 dB SPL
- - - Magnetic input: 31.6 mA/m

Thin tube 0.9
— Acoustic input: 60 dB SPL
- - - Magnetic input: 31.6 mA/m


	Hook	Thin tube 1.3 mm	Thin tube 0.9 mm
OSPL90, Peak (dB SPL)	123	119	117
OSPL90, 1600 Hz (dB SPL)	121	114	108
OSPL90, HFA (dB SPL)	119	115	110
Full-on gain, Peak (dB) ¹	54	54	56
Full-on gain, 1600 Hz (dB) ¹	48	48	44
Full-on gain, HFA (dB) ¹	48	48	44
Reference test gain (dB)	42	37	34
Frequency range (Hz)	100-7300	100-6300	100-6800
Telecoil output, 1 mA/m field (1000 Hz) (dB SPL)	75	85	84
Telecoil output, HFA SPLITS L/R (dB SPL)	100	97	91
Total harmonic distortion (Input 70 dB SPL), 500 Hz (%)	<4	<4	<2
Total harmonic distortion (Input 70 dB SPL), 800 Hz (%)	<3	<2	<2
Total harmonic distortion (Input 65 dB SPL), 1600 Hz (%)	<2	<2	<2
Equivalent input noise level, Omni (dB SPL)	17	19	21
Equivalent input noise level, Dir (dB SPL)	29	*	*
Battery consumption, Typical (mA) ²	2.0	1.9	2.0
Battery consumption, Quiescent (mA) ²	1.9	1.9	1.9
Battery life, artificial measurement, hours ³	90	95	90
Expected battery life, hours (battery size 312 – IEC PR41) ⁴	50-55	*	*

* No measurement performed


- Measured with the gain control of the hearing aids set to their full-on position minus 20 dB and with an input SPL of 70 dB. This is to obtain a gain response equal to the full-on gain response from e.g. IEC 60118-0:1983+A1:1994 but without influence of feedback.
- Battery current is measured according to IEC 60118-0:1983/AMD1:1994 §7.11, IEC 60118-0:2015 §7.7 and ANSI S3.22:2014 §6.13 after a settling time of minimum 3 minutes.
- Based on the standardised battery consumption measurement (IEC 60118-0:1983/AMD1:1994). The actual battery life depends on battery quality, use pattern, active feature set, hearing loss and sound environment.
- Real usage battery life is shown as an estimated interval based on mixed use cases with variable amplification settings and variable input levels, incl. direct stereo streaming from a TV (25% of the time) and streaming from a mobile phone (6% of the time).

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hearingsolutions.philips.com

Philips and the Philips Shield Emblem are registered trademarks of Koninklijke Philips N.V. and are used under license. This product has been manufactured by or for and is sold under the responsibility of SBO Hearing A/S, and SBO Hearing A/S is the warrantor in relation to this product.