TEKNİK VERİLER

Entra B 20 | 10

ITC, ITE HS, ITE FS

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 SBO Hearing A/S
Kongebakken 9
DK-2765 Smørum
Danimarka

www.bernafon.com

Bernaфон is part of the Demant Group.

TECHNICAL DATA

Entra B 20 | 10

ITC, ITE HS, ITE FS

Bernafon Entra B ITC, ITE HS and ITE FS are Bernafon's most flexible in-the-ear hearing instruments, suitable for slight to severe hearing losses. These instruments offer the revolutionary Hybrid Technology™ included in other Bernafon

hearing instrument styles. Featuring NFMI and optional Bluetooth® Low Energy technology, they can be used to stream audio directly to the hearing instruments.

Speaker 75 | 90 | 100

ITC (In-The-Canal)

Speaker 75 | 90 | 100

ITE HS (Half Shell)

Speaker 75 | 90 | 100

ITE FS (Full Shell)

Technical features

- Hands-free communication^{1,2}
- Direct audio streaming^{1,3}
- Directional microphones
- NFMI (Near-Field Magnetic Induction)
- Bluetooth Low Energy technology (optional)
- Push-button (optional)
- Volume control (optional)
- Telecoil (optional)

Accessories¹

- Bernafon App
- RC-A (remote control)
- TV-A (TV adapter)
- SoundClip-A

General features

- Digital Programmable
- Automatic or Manual Volume Control
- Maximum Output Control System
- MPO-Maximum Power Output
- GC-Gain Control
- AGC-Automatic Gain Control
- Noise Reduction
- Feedback Management
- Dual Microphone
- 4 Programs (optional)
- FM Compatible(with Telecoil)
- 48 Channels

For information on compatibility, please visit www.bernafon.com/compatibility

Operating and charging conditions

Temperature: +1°C to +40°C (34°F to 104°F)
Humidity: 5% to 93% relative humidity, non-condensing
Atmospheric pressure: 700 hPa to 1060 hPa

Transportation and storage conditions

Temperature and humidity shall not exceed the mentioned limits for extended periods during transportation and storage:

Transport

Temperature: -25°C to + 60°C (-13°F to 140°F)
Humidity: 5% to 93% relative humidity, non-condensing
Atmospheric pressure: 700 hPa to 1060 hPa

Storage

Temperature: -25°C to + 60°C (-13°F to 140°F)
Humidity: 5% to 93% relative humidity, non-condensing
Atmospheric pressure: 700 hPa to 1060 hPa

1) Only available for hearing aids with Bluetooth Low Energy technology

2) Hands-free communication is available on select devices

3) From iPhone, iPad, Mac and select Android devices

WARNING: No modification of this equipment is allowed.

Apple, the Apple logo, iPhone, iPad, Mac and the Mac logo are trademarks of Apple Inc., registered in the U.S. and other countries. Use of the Made for Apple badge means that an accessory has been designed to connect specifically to the Apple product(s) identified in the badge, and has been certified by the developer to meet Apple performance standards. Apple is not responsible for the operation of this device or its compliance with safety and regulatory standards. Android™ is a trademark of Google LLC.

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


Feature overview

	Entra B 20	Entra B 10
Hybrid Technology™		
Hybrid Sound Processing™	•	•
Frequency bandwidth	8 kHz	8 kHz
Speech Balancer	•	•
Hybrid Noise Management™	•	•
Smart Noise Reduction	3 options	2 options
Smart Directionality	4 options	3 options
Hybrid Feedback Canceller™	•	•
Speech		
Low Frequency Enhancer ¹	•	•
Frequency Composition™	•	•
Comfort		
Transient Noise Reduction	2 options	—
Wind Noise Manager	•	•
Soft Noise Manager	•	•
Directionality controls		
Dynamic	•	—
Adaptive Full Directionality	•	•
Fixed Directionality	•	•
Fixed Omni	•	•
Individualisation		
Personalisation	•	•
Fitting bands	14	12
Program options ² / memories ³	10/4	8/4
Music Experience	•	—
Binaural coordination: VC, program change ³	•	•
Automatic Adaptation Manager	•	•
Transition	•	•
Data Logging	•	•
Tinnitus SoundSupport ⁴	•	•


1) Requires Bluetooth Low Energy technology
2) Can vary if no telecoil present
3) Requires either Bluetooth Low Energy technology or push-button
4) Requires push-button

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

Technical information
Omnidirectional mode is used
unless otherwise stated.

Speaker 75 / 90
— Acoustic input: 60 dB SPL
- - - Magnetic input: 31.6 mA/m

Speaker 100
— Acoustic input: 60 dB SPL
- - - Magnetic input: 31.6 mA/m

OSPL90, Peak (dB SPL)	120
OSPL90, 1600 Hz (dB SPL)	110
OSPL90, HFA (dB SPL)	111
Full-on gain, Peak (dB) ¹	58
Full-on gain, 1600 Hz (dB) ¹	48
Full-on gain, HFA (dB) ¹	48
Reference test gain (dB)	36
Frequency range (Hz)	100-7500
Telecoil output, 1 mA/m field (1600 Hz) (dB SPL)	79
Telecoil output, 10 mA/m field (1600 Hz) (dB SPL)	99
Total harmonic distortion (Input 70 dB SPL), 500 Hz (%)	<2
Total harmonic distortion (Input 70 dB SPL), 800 Hz (%)	<2
Total harmonic distortion (Input 70 dB SPL), 1600 Hz (%)	<3
Equivalent input noise level, Omni (dB SPL)	18
Equivalent input noise level, Dir (dB SPL)	26
Battery consumption, Typical (mA) ²	1.9
Battery consumption, Quiescent (mA) ²	1.9
Battery life, artificial measurement, hours ³	95
Expected battery life, hours (battery size 312 – IEC PR41) ⁴	55-60

Speaker 75

Speaker 90

Speaker 100

OSPL90, Peak (dB SPL)	129	134
OSPL90, 1600 Hz (dB SPL)	124	134
OSPL90, HFA (dB SPL)	124	131
Full-on gain, Peak (dB) ¹	63	72
Full-on gain, 1600 Hz (dB) ¹	60	70
Full-on gain, HFA (dB) ¹	59	67
Reference test gain (dB)	49	60
Frequency range (Hz)	100-7500	100-7500
Telecoil output, 1 mA/m field (1600 Hz) (dB SPL)	90	101
Telecoil output, 10 mA/m field (1600 Hz) (dB SPL)	110	121
Total harmonic distortion (Input 70 dB SPL), 500 Hz (%)	<2	<2
Total harmonic distortion (Input 70 dB SPL), 800 Hz (%)	<3	<3
Total harmonic distortion (Input 70 dB SPL), 1600 Hz (%)	<2	<3
Equivalent input noise level, Omni (dB SPL)	15	12
Equivalent input noise level, Dir (dB SPL)	24	23
Battery consumption, Typical (mA) ²	2.1	2.0
Battery consumption, Quiescent (mA) ²	1.9	1.9
Battery life, artificial measurement, hours ³	85	90
Expected battery life, hours (battery size 312 – IEC PR41) ⁴	40-60	50-60


- Full-on gain is measured with the gain control of the hearing instruments set to its full-on position minus 20 dB and with an input SPL of 70 dB. This is to obtain a gain response equal to the full-on gain response from e.g. IEC 60118-0+A1:1994 but without influence of feedback.
- Battery current is measured according to IEC 60118-0:1983/AMD1:1994 §7.11, IEC 60118-0:2015 §7.7 and ANSI S3.22:2014 §6.13 after a settling time of minimum 3 minutes.
- Based on the standardised battery consumption measurement (e.g. IEC 60118-0:1983/AMD1:1994). The actual battery life depends on battery quality, use pattern, active feature set, hearing loss and sound environment.
- Real usage battery life is shown as an estimated interval based on mixed use cases with variable amplification settings and variable input levels, not direct stereo streaming from a TV (25% of the time) and streaming from a mobile phone (6% of the time).

İSBU BELGE (ASLINDAN, FOTOKOPİ - FAX, İMZA) NİNDEN
TARAFIMIZIN İNGİLİZCE DEN TÜRKÇE'YE ÇEVİRİLMİŞTİR.
Emrah B. YEMİNLİ MÜTERCİM

Entra B 20 | 10 ITC/ITE HS/FS

2CC Coupler



Measured according ANSI S3.22: 2014,
IEC 60118-0:2015 and IEC 60318-5:2006.

Technical information
Omnidirectional mode is used
unless otherwise stated.

Speaker 75 / 90
— Acoustic input: 60 dB SPL
- - - Magnetic input: 31.6 mA/m

Speaker 100
— Acoustic input: 60 dB SPL
- - - Magnetic input: 31.6 mA/m

	Speaker 75	Speaker 90	Speaker 100
OSPL90, Peak (dB SPL)	109	119	126
OSPL90, 1600 Hz (dB SPL)	102	115	126
OSPL90, HFA (dB SPL)	103	116	123
Full-on gain, Peak (dB) ¹	46	54	64
Full-on gain, 1600 Hz (dB) ¹	40	51	63
Full-on gain, HFA (dB) ¹	40	51	60
Reference test gain (dB)	26	39	46
Frequency range (Hz)	100-7500	100-7500	100-5400
Telecoil output, 1 mA/m field (1000 Hz) (dB SPL)	69	80	91
Telecoil output, HFA SPLITS L/R (dB SPL)	85	98	105
Total harmonic distortion (Input 70 dB SPL), 500 Hz (%)	<2	<2	<2
Total harmonic distortion (Input 70 dB SPL), 800 Hz (%)	<2	<2	<2
Total harmonic distortion (Input 65 dB SPL), 1600 Hz (%)	<2	<2	<2
Equivalent input noise level, Omni (dB SPL)	17	15	15
Equivalent input noise level, Dir (dB SPL)	27	27	30
Battery consumption, Typical (mA) ²	2.0	2.4	2.1
Battery consumption, Quiescent (mA) ²	1.9	1.9	1.9
Battery life, artificial measurement, hours ³	90	75	85
Expected battery life, hours (battery size 312 – IEC PR41) ⁴	55-60	40-60	50-60

- Full-on gain is measured with the gain control of the hearing instruments set to its full-on position minus 20 dB and with an input SPL of 70 dB. This is to obtain a gain response equal to the full-on gain response from e.g. IEC 60118-0+A1:1994 but without influence of feedback.
- Battery current is measured according to IEC 60118-0:1983/AMD1:1994 §7.11, IEC 60118-0:2015 §7.7 and ANSI S3.22:2014 §6.13 after a settling time of minimum 3 minutes.
- Based on the standardised battery consumption measurement (e.g. IEC 60118-0:1983/AMD1:1994). The actual battery life depends on battery quality, use pattern, active feature set, hearing loss and sound environment.
- Real usage battery life is shown as an estimated interval based on mixed use cases with variable amplification settings and variable input levels, incl. direct stereo streaming from a TV (25% of the time) and streaming from a mobile phone (6% of the time).

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 SBO Hearing A/S
Kongebakken 9
DK-2765 Smørum
Danimarka

www.bernafon.com

Bernafon is part of the Demant Group.

TECHNICAL DATA

Entra B 20 | 10

IIC, CIC

Bernafon Entra B IIC and CIC are the smallest in-the-ear hearing instruments of the Entra B family, suitable for slight to severe hearing losses. These instruments offer the revolutionary Hybrid Technology™ included in other Bernafon instrument styles. Sophisticated features work together for

seamless and boundless adaptation to listening environments. Placed deeply in the canal, these instruments allow the user to benefit from the natural pinna effect and improve their sound localisation abilities. Each style supports two power levels to better accommodate users' needs.

Technical features

- Battery size: 10
- Hydrophobic coating
- Push-button¹
- NFMI (Near-Field Magnetic Induction)¹

General features

- Digital Programmable
- Automatic Volume Control
- Maximum Output Control System
- MPO-Maximum Power Output
- GC-Gain Control
- AGC-Automatic Gain Control
- Noise Reduction
- Feedback Management
- Single Microphone
- 4 Programs (only optional for CIC)
- 48 Channels

Operating and charging conditions
Temperature: +1°C to +40°C (34°F to 104°F)
Humidity: 5% to 93% relative humidity,
non-condensing
Atmospheric pressure: 700 hPa to 1060 hPa

Transportation and storage conditions
Temperature and humidity shall not exceed the
mentioned limits for extended periods during
transportation and storage.

Transport
Temperature: -25°C to +60°C (-13°F to 140°F)
Humidity: 5% to 93% relative humidity,
non-condensing
Atmospheric pressure: 700 hPa to 1060 hPa

Storage
Temperature: -25°C to +60°C (-13°F to 140°F)
Humidity: 5% to 93% relative humidity,
non-condensing
Atmospheric pressure: 700 hPa to 1060 hPa

1) Optional features only available for CIC

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

1) Requires NFMI (only available for CIC)

2) Requires push-button (only available for CIC)

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





- 1) Full-on gain is measured with the gain control of the hearing instruments set to its full-on position minus 20 dB and with an input SPL of 70 dB. This is to obtain a gain response equal to the full-on gain response from e.g. IEC 60118-0+A1:1994 but without influence of feedback.
- 2) Battery current is measured according to IEC 60118-0:1983/AMD1:1994 §7.11, IEC 60118-0:2015 §7.7 and ANSI S3.22:2014 §6.13 after a settling time of minimum 3 minutes.
- 3) Based on the standardised battery consumption measurement (e.g. IEC 60118-0:1983/AMD1:1994). The actual battery life depends on battery quality, use pattern, active feature set, hearing loss and sound environment.
- 4) Real usage battery life is shown as an estimated interval based on mixed use cases with variable amplification settings and variable input levels.

Entra B 20 | 10 IIC / CIC

2CC Coupler

Measured according to ANSI S3.22:2014,
IEC 60118-0:2015 and IEC 60318-5:2006

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